

SAMPLE SUBMISSION GUIDE

Human Cell Line Authentication Service by STR/DNA Profiling – 2.0

Service

Cell Line Authentication (CLA) Service – Basic Service (from genomic DNA) – incl. Prepaid Option

- Cell Line Authentication Service for human cell lines is performed using Applied Biosystems™ AmpFLSTR™ Identifier™ Plus PCR Amplification Kit with 16 markers (CSF1PO, D2S1338, D3S1358, D5S818, D7S820, D8S1179, D13S317, D16S539, D18S51, D19S433, D21S11, FGA, TH01, TPOX, vWA, and the gender marker Amelogenin) for high-resolution discrimination and includes all relevant markers of ANSI/ATCC ASN-0002-2011. Authentication of human cell lines: Standardization of STR profiling.

Cell Line Authentication (CLA) Service – Additional Services

- DNA Extraction from cell pellets.
- Comparison of the DNA profile (database entries e.g. at DSMZ, ATCC, etc.*) plus official certificate (accepted by peer-reviewed journals).
- Chromatogram of the DNA profile as PDF.

Results: Certificate will be uploaded into your account in a turnaround time (TAT) of 5 working days.

*Is the reference genetic fingerprint of your cell line just available in a special database or in a publication, please provide us with this information directly at the ordering step.

Shipment Instructions

- Order CLA Prepaid Barcodes for the Basic Service.
- Click on the '[New Order Page](#)' Button on the CLA page or on '[Cell Line Authentication 2.0](#)' in our Order Menu.
- Follow the steps to register your samples and prepaid barcodes.
- Send your samples directly to below address or place them into one of our dropboxes (+approx. one day logistics time).
- Samples are accepted as:
 - Cell Pellets (preferred) 2 mL reaction tube (no falcon tubes) with >10.000 cells (**visible cell pellet!**)
Cell culture medium completely removed; cell pellets washed with TE buffer or PBS; shipped at ambient temperature or cooled (no need for dry ice)
 - Genomic DNA 2 mL reaction tube; volume >20 µL; amount >10 ng/µl
shipped at ambient temperature
- If >20 samples per batch are shipped, please sort the 2 mL tubes in a box.

- **Shipment Address:** **APG Department**
Eurofins Genomics
Anzinger Str. 7a
85560 Ebersberg, Germany

Contact: Please do not hesitate to contact us in case of any questions!

Email: ngs-support@eurofins.com

Phone: +49 7531 81 60 68 (Support)

Registered trademarks listed at <http://www.eurofinsgenomics.eu>. The laboratory is accredited acc. to DIN EN ISO/IEC 17025:2005. The accreditation applies only to the test methods specified in the accreditation certificate. All analyses have been carried out with greatest care and on the basis of state of the art scientific knowledge. The results refer solely to the analysed samples. The duplication and publication also in parts requires a written authorization by this laboratory. Our General Terms and Conditions apply exclusively and are available under eurofinsgenomics.com.

