

NGS Barcodes

Sample Submission Information for NGS Products

IMPORTANT UPDATE FOR NGS SAMPLE SHIPMENTS

Order first, then send samples –

Samples without a corresponding order cannot be registered and may, in the worst case, be discarded. Therefore, always place your order first before sending samples.

Sample requirements


Sample volume, concentration, and input depend on the product. Always check the Sample Submission Guide for the latest requirements. NanoDrop values may be inaccurate; Qubit is preferred.

IMPORTANT SHIPPING INFORMATION

- Please use 1.5 ml snap-cap tubes (safe-lock)
- Stick barcodes on the outside of the tube and in correct order (Horizontal not vertical)
- For projects with > 20 samples, please sort the samples by Barcode and number them on the lid
- Samples for extraction need to be shipped with the risk group classification sheet
- Frozen barcodes tend to fall apart. Please ship in cryoboxes or other material that keeps barcodes close to the sample. Number samples on the lid
- If you use dry ice please store samples in a separate bag or box

Please send your samples to:

Eurofins Genomics FEB
Anzinger Str. 7a
85560 Ebersberg / Germany



Sample Preparation & Submission.

» **Frozen barcodes** tend to fall apart. Please ship in cryoboxes or other material that keeps barcodes close to the sample. Number samples on the lid.

1.




1. DEFINE NGS PROJECT

For direct orders please follow the instructions on the respective order page.

For quotes please click on “Complete sample information” in my quotes and fill in the required details.

2.



2. LABEL TUBES

Affix the barcode on your tube as shown here.

Please use 1.5 ml snap-cap tubes (safe-lock).

Do not use parafilm to seal your tubes.

3.



3. SUBMIT SAMPLES

Please add the order confirmation and/or your sample submission form to your order.

4.



4. STATUS TRACKING & RESULTS

The status of your order can be tracked online.

You will receive an Email once your data is ready to download.

Additional Service & Information

» Please read our [sample submission guide](#) carefully

» Looking for a prepaid option for NGS? Learn more about our [NGS Coupons](#)

» Unsure about the best product choice? Send us your request

» Need UPS labels for shipping? [Request here](#)



A COMPREHENSIVE OVERVIEW OF OUR NGS PORTFOLIO, INCLUDING OUR OXFORD NANOPORE SERVICES, CAN BE FOUND ON OUR [WEBSITE](#)

